

## New Product

## Syringe cap for enteral feeding syringes



The multifunctional LOGICA Syringe cap is used to hygienically seal enteral feeding syringes with a filling volume of 20 ml or 60 ml.

Due to the combined design, the cap can be used either to close the 20 ml or the 60 ml syringe barrel. It is no longer necessary to stock two different sizes.

The built-in ventilation filter enables hygienic ventilation of the feeding syringes during gravity application.

- > hydrophobic filter material
- > latex-free
- > PVC-free
- > DEHP-free

---

**Closure for 20 ml syringe barrel**

compatible to 20 ml enteral feeding syringes  
 (e.g. Fresenius® Freka® Connect ENFit™ / ProNeo 20 ml)

---

**Closure for 60 ml syringe barrel**

compatible to 60 ml enteral feeding syringes  
 (e.g. Fresenius® Freka® Connect ENFit™ / ProNeo 60 ml)

---

**Hydrophobicity of filter**

> 93 mmHg (> 124 mbar)

---

**Pore size of filter**

5 µm

---

**Sterilization Compatibility**

ETO-gas,  $\gamma$ -irradiation and  $\beta$ -irradiation up to 50 kGy (cumulative)

---

<sup>\*)</sup>: Closure diameters can be adapted to customer-specific syringe barrels.